


HARNESS SOLUTIONS

CAREVA HARNESS SYSTEM



CAREVA COMBINATION HARNESS




HARNESS SOLUTIONS

CAREVA HARNESS SYSTEM



CAREVA COMBINATION HARNESS




The **Careva Harness System** is a postural and support harness that provides upper body support and prevents the person falling forward or to the side to ensure that the vehicles seatbelt works correctly.

What is the function of **Careva Combination Harness**?

This product is a full-body support, that is used to ensure the occupant is safely seated in the correct seating position for the vehicles seatbelt.

The system consists of a range of single belts that can be combined as required for users of all ages who need additional support for an upright sitting position. **Careva Combination Harness** can be used in any seat, provided there is an opening between the seat and the backrest.



FEATURES

- Universal postural belt system that can be used in any seat, provided there is an opening between the seat and the backrest
- Quick and easy to install and transfer between vehicles

SPECIFICATIONS

- Available in: **Small** 4-7 years
Medium 8-12 years
Large 12 years - Adult
- Can be used with wheelchairs with a stable backrest
- Careva Harness System DOES NOT** replace the vehicle safety belt

REQUIREMENTS It is a requirement to keep a **current** Doctor's Certificate in your vehicle when this Harness is in use in the vehicle

The **Careva Combination Harness** is a positioning belt designed with the help of many involved parents of children with disabilities, drivers in the community transport, Occupational Therapists and health professionals to ensure the correct postural supports are achieved.

Our harnesses are quick to install and can be easily transferred from across vehicles.

BENEFITS


- Ensures the occupants safety while travelling in the vehicle
- No modifications to your vehicle
- Provides you with transport independence

OPTIONAL ACCESSORIES



THORACIC BELT

An additional option that provides additional upper body chest support



4/45 Salisbury Road
ASQUITH NSW 2077
(02) 9482 4572



www.mobilityengineering.com.au

sales@mobilityengineering.com.au