

EASY TRAVEL MAGNET HARNESS

The Easy Travel Magnet Harness is a behavioural support harness.

The Easy Travel Magnet Harness is used to ensure the occupant is in the correct seating position for the vehicles seatbelt to work correctly. This product is a vest designed with magnet buttons down the front of the harness, which cannot be undone without the magnetic key. The design of this harness is for people with behaviours of concern and have the tendency to get out of the vehicles seating position.

The Easy Travel Magnet Harness is to be used in accordance with the vehicle's original seat and seatbelt. It can be installed into a variety of vehicles without any modifications and can be transferred across different vehicles. This harness can also be worn before getting into the vehicle if required, and can be installed onto seats using either vertical or horizontal straps around the seatback.



	Measurement A	Measurement B
Extra Small	37cm	45-60cm
Small	42cm	55-70cm
Medium	50cm	60-85cm
Large	55cm	65-95cm
Extra Large	59cm	70-100cm

Product Features:

- Universal behavioural and postural belt system that can be used in any seat
- Quick and easy to install and transfer between vehicles
- Can only be opened with the included magnet key, or a conventional magnet
- Machine washable at 40°C
- Comfortable design for the occupant
- No modifications to the vehicle required for use.

Quick Measurement Guide

XS-XL sizes available. See above table.

Important Note's

- Harness **does not** replace the vehicle safety belt
- Anti-submarining strap as standard
- It is a requirement to keep a **current** doctor's certificate in your vehicle when this harness is in use in the vehicle.