

# HARNESS SOLUTIONS

## B&S 7-POINT HARNESS



# HARNESS SOLUTIONS

## B&S 7-POINT HARNESS






This system is used for children and young adults who often try to remove themselves from their harness system. The harness consists of a baseband that is easy to attach to the rear seat of a vehicle and includes a fully adjustable harness with 7 lockable buckles. As an added bonus, the strap between the strap can be removed when not in use.



In the event that multiple vehicles are used, an additional, separately supplied base band can be incorporated. The universal harness is used in both vehicles and the base belts can remain attached.

The B&S Car restraint system can be used in combination with a booster seat. This must be a booster seat without a backrest. It is, of course, important that you keep the restraint vest at the correct height. An extra point of attention is the 7th point, the belt that goes between the legs. This must pass behind the booster seat!




### FEATURES

-  Can be used in conjunction with booster seats.
-  The leg strap can be optional.
-  Safety lock with release pin.

### SPECIFICATIONS

-  Complies with recommended safety standards.
-  Recommend for children and young adults.

### BENEFITS

-  Safely secures the occupant.
-  Can be used across multiple vehicles or buses.
-  Provides a total of 7 safety straps for security.

